## CH1202: Lab Report II

# Determination of Degree of Hydrolysis and Hydrolysis Constant by Potentiometry

## ABHISRUTA MAITY

#### 21MS006

am21ms006@iiserkol.ac.in

Indian Institute of Science Education and Research, Kolkata Mohanpur, West Bengal, 741246, India

#### **Contents**

1	Aim	1
2	Apparatus Required	1
3	Chemicals Required	1
4	Experimental Data	2
5	Conclusion	2

# §1 Aim

To determine the degree of hydrolysis and hydrolysis constant of  $Anilinium\ Hydrochloride$  using Potentiometer.

# §2 Apparatus Required

- Potentiometer
- Platinum Electrode
- Calomel Electrode

# §3 Chemicals Required

- Anilinium Hydrochloride
- Quinhydrone

Lab Report II 5 Conclusion

### §4 Experimental Data

#### Calculation of Hydrolysis Constant

1. pH is given by

$$pH = \frac{-E_{\text{obs}} + E_{\text{QH}} + E_{\text{cal}}}{0.0591}$$

where  $E_{QH} = 0.6996 V$  and  $E_{cal} = -0.242 V.^{1}$ 

- 2. Since  $pH = -\log(H^+) = -\log(c\alpha)$ ,  $pH = -\log(c) \log(\alpha)$ , the degree of hydrolysis  $\alpha$  can be calculated at any given concentration.
- 3. From  $\alpha$ , using  $K_h = \frac{c\alpha^2}{1-\alpha}$ , we can deduce hydrolysis constant of Anilinium Hydrochloride.
- 4. The dissociation constant can also be calculated using the relation  $K_b = \frac{K_w}{K_h}$ .

$\left[\mathrm{C_6H_5NH^{3+}Cl^{-}}\right]$	$E_{\rm obs}$	pH	$\alpha \ [\times 10^{-2}]$	$K_h \ [\times 10^{-5}]$	$K_b \ [\times 10^{-10}]$
0.10	0.277	3.06	0.88	0.78	12.8
0.05	0.274	3.11	1.56	1.24	8.04
0.02	0.264	3.28	2.65	1.44	6.93
0.01	0.260	3.34	4.53	2.15	4.64

## §5 Conclusion

The experimented value of  $K_h = 1.40 \times 10^{-5}$  and that of  $K_b = 8.11 \times 10^{-10}$ .

 $<sup>^1{\</sup>rm These}$  are oxidation potentials. We took the sign conventions accordingly.

<sup>&</sup>lt;sup>2</sup>The value of  $K_w$  at 25° C is assumed as 10<sup>-14</sup>.